



Vectors

[In vivo Biotinylation-Strains](#)

[In vitro Biotinylation-Enzyme](#)

[Purification & Detection](#)

[Biotinylated recombinant proteins](#)



[Mammalian Cell-Free Translation System](#)

[Protocols](#)

[Literature](#)

[FAQs](#)

[New Products](#)

Protocols

Avidity provides protocols for guided use of products utilized in AviTag™ Technology. If further assistance is needed, please contact our technical support at **1-877-333-6063**.

[Induction of bacterial strains for in vivo biotinylation of AviTag™ proteins](#)

[General reaction conditions for biotin ligase](#)

[In vitro labeling purified AviTag™ proteins with biotin ligase](#)

[In vitro labeling AviTag™ proteins in BACTERIAL EXTRACTS](#)

[In vitro labeling AviTag™ proteins in MAMMALIAN EXTRACTS](#)

[In vitro labeling AviTag™ proteins on CELL SURFACES](#)

[In vitro labeling MHC MOLECULES with AviTag™](#)

[Guidelines for Anti C-terminus AviTag monoclonal antibody \(AbCA\)](#)

[Determine the extent of biotinylation with BIS-300 \(BRTA assay\)](#)

[Standard Protocol for AvidExpress™ cell-free translation system](#)